



PHARMACY NETWORK SERVICES

New Press-In Bottle Adapters for Lorazepam Solution

Lorazepam Solution leaking on you? Is this creating a mess in the carts and on the controlled substance disposition sheet? Well PNS has come up with a solution! New Press-In Bottle Adapters. You may have used these before, but this should solve the problem with lorazepam solution leaking!

What does it look like?



The pharmacy will have already inserted these prior to receiving the lorazepam solution. Here is what it will look like with the bottle:



Be sure to insert the syringe prior to inverting the bottle to draw back.

Please let us know if you have any issues or concerns, we are here to help!